

## **Persyst 15Rev.C1 (build 2025.12.09) Marketing Release Notes**

The latest update from Persyst focuses on improving **workflow reliability, stability, and data integrity across the Persyst ecosystem and key partner integrations.**

### Persyst Trending

- Resolved an issue causing Trend files to grow excessively and exceed file size limits when a second record was accessed during a live recording.

### Persyst EEG Viewer

- Resolved a crash that occurred when closing the application.

### Persyst ESI

- Removed an unnecessary step during MRI selection that required user intervention and prevented the upload progress bar from starting immediately.

### Nihon Kohden

- Resolved a file path handling issue that prevented files from loading correctly when launched from the NK Polaris interface.

### Alpha Trace

- The Alpha Trace data reader was updated to correctly import customer channel names and respect channel activation flags.